

Product Name :
Bedside Locker Deluxe

Product Code :
NLABMEDICAL-M2500005



Description :

Bedside Locker Deluxe

Technical Specification :

Specification :-

Frame Work of Bedside Locker Deluxe is made of mild steel pressed sheet around square tubes.

It has stainless steel top with raise edges.

It is closed from three sides.

It comes with cupboard box and one drawer for storage.

Bedside Locker Deluxe is Pre Treated Epoxy Powder Coated.

Product Size : 40 L x 40 W x 82 H cms.

We are Manufacturer, Exporters and Suppliers of Bedside Locker Deluxe at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003