

Product Name :
Bedside Locker Laminated Wood

Product Code :
NLABMEDICAL-M2500004



Description :

Bedside Locker Laminated Wood

Technical Specification :

Specification :-

Frame Work of Bedside Locker with Cabinet is made of mild steel pressed sheet.

It has Wooden Top, Wooden door and Wooden drawer face.

It comes with cupboard box and one drawer for storage.

Bedside Locker with cabinet is Pre Treated Epoxy Powder Coated.

Product Size : 40 L x 40 W x 82 H cms.

We are Manufacturer, Exporters and Suppliers of Bedside Locker Laminated Wood at tender and export competitive prices from India.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments