

**Product Name :**  
Bedside Locker Standard Plain Top

**Product Code :**  
NLABMEDICAL-M2500001



**Description :**

Bedside Locker Standard Plain Top

**Technical Specification :**

Frame Work of Bedside Locker Standard is made of a Square tube.

It has a Stainless Steel plain top.

It comes with a cupboard box for storage.

Bedside Locker Standard is Pre Treated Epoxy Powder Coated.

Product Size : 40 L x40 W x 82 H cms.

We are Manufacturers, Exporters, and Suppliers of Bedside Locker Standard Plain Top at tender and export competitive prices from India.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

