

**Product Name :**  
Bedside Locker Standard Tray Top

**Product Code :**  
NLABMEDICAL-M2500002



**Description :**

Bedside Locker Standard Tray Top

**Technical Specification :**

Specification :-

Frame Work of Bedside Locker Standard is made of Square tube.

It has Stainless Steel tray top.

It comes with cupboard box for storage.

Bedside Locker Standard is Pre Treated Epoxy Powder Coated.

Product Size : 40 L x40 W x 82 H cms.

We are Manufacturer, Exporters and Suppliers of Bedside Locker Standard Tray Top at tender and export competitive prices from India.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments