

**Product Name :**  
Bedside Locker with Overbed Table

**Product Code :**  
NLABMEDICAL-M2500007



**Description :**

Bedside Locker with Overbed Table

**Technical Specification :**

Specification :-

Frame Work of Bedside Locker is made of mild steel pressed sheet around square tubes.

It has laminated cabinet door, drawer face and top.

The overbed table is laminated with height adjustable from 90 upto 115cm.

It has bottle holder at bottom.

Product Size : 40 L x 40 W x 90 H cms.

We are Manufacturer, Exporters and Suppliers of Bedside Locker with Overbed Table at tender and export competitive prices from India.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments