

Product Name :
Bein Root Elevator

Product Code :
NLABMEDICAL-M2340001



Description :

Bein Root Elevator

Technical Specification :

Bein Root Dental Elevator is made of High-Grade Stainless Steel.

Dental Elevators are used to luxate teeth and this technique imparts forces to tooth particles that sever the periodontal ligament around tooth roots inside the socket and expand alveolar bone around tooth particles. These effects can result in extraction of the tooth particles.

We are Manufacturer, Exporters and Suppliers of Bein Root Elevator at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003