

Product Name :
Bell Crank Lever

Product Code :
NLAB-MECHANICALAB240032



Description :

Bell Crank Lever

Technical Specification :

A wooden right angled lever with fulcrum at the angle, the horizontal arm is of 15 cm. and vertical arm is of 20 cm., a spring balance with adjustable wing nut is attached to the shorter arm. The lever arm is with scale and grooved at every 5 cm. Complete with one sliding weight of 1 kg.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003