

Product Name :
Bell in Vacuum

Product Code :
NLAB-EDUSCHOOLEQ155002



Description :

Bell in Vacuum

Technical Specification :

A versatile unit for demonstrating the properties of a vacuum.

Use with hand pump or ideally with the electric vacuum pump.

Comprises a sturdy baseplate with the tap, glass bell jar and detachable battery operated bell

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003