

Product Name :
Belt Flex Tester Machines

Product Code :
NLAB-TESTINGLABM98001



Description :

Belt Flex Tester Machines

Technical Specification :

Belt Flex Tester. We use high end quality components and metals to make these equipments following industry set guidelines. Owing to their outstanding performance and high efficiency, Belt Flex Tester is used to flex the soles of footwear with their actual patterns

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003