

**Product Name :**  
Belt Sander

**Product Code :**  
NLAB-MACHINELABWK33016



**Description :**

Belt Sander

**Technical Specification :**

Oscillating

Power: 1470W

Belt speed: 960m/min

Belt size: 150 x 2260mm

Table size: 680 x 305mm

Dust extraction outlet: 2 x 100mm

Wood Workshop Portable Tools Belt Sander Manufacturer, Supplier and Belt Sander Exporter in India

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003