

Product Name :
Bench Model Heavy Duty Tapping Machine 5MM - 12MM

Product Code :
NLAB-MACHINELABWK12007



Description :

Bench Model Heavy Duty Tapping Machine 5MM - 12MM

Technical Specification :

Bench Model Heavy Duty Tapping Machine 5MM - 12MM .

PARAMETRICULARS

- T**apping capacity in Steel
- T**apping capacity in Brass & Aluminium
- T**able 145
- D**iameter of column
- S**pindle Centre to column
- D**istance chuck to work table
- T**apping Depth
- S**ize of base plate L x W
- T**hroat (Z Saction)
- S**peed of spindle speeds

1400 Rev. R.P.M. H.P.
25 Fan 55
121 Weight Approx.
676g. Weight Approx.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003