

**Product Name :**  
Bench Model Heavy Duty Tapping Machine 5MM - 12MM

**Product Code :**  
NLAB-MACHINELABWK12007



**Description :**

Bench Model Heavy Duty Tapping Machine 5MM - 12MM

**Technical Specification :**

Bench Model Heavy Duty Tapping Machine 5MM - 12MM .

**PARAMETRICULARS**

- T**apping capacity in Steel
- T**apping capacity in Brass & Aluminium
- T**able 145
- D**iameter of column
- S**pindle Centre to column
- D**istance chuck to work table
- T**apping Depth
- S**ize of base plate L x W
- T**hroat (Z Saction)
- S**peeds of spindle speeds

1400 Rev. R.P.M. H.P.  
25 Fan 55  
Weight Approx.  
Weight Approx.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003