

Product Name :
Bench Top Mechanical Heat Pump Trainer**Product Code :**
REFI0003**Description :**

Bench Top Mechanical Heat Pump Trainer

Technical Specification :

Bench Top Mechanical Heat Pump Trainer: refrigeration and air conditioning lab equipment: air-to-water heat pump trainer consists of a compressor, an evaporator with fan, a thermostatic expansion valve and a coaxial coil heat exchanger as condenser.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773