

**Product Name :**  
Bench Top Mechanical Heat Pump Trainer

**Product Code :**  
NLAB-MECHANICALAB170003



**Description :**

Bench Top Mechanical Heat Pump Trainer

**Technical Specification :**

Bench Top Mechanical Heat Pump Trainer: refrigeration and air conditioning lab equipment: air-to-water heat pump trainer consists of a compressor, an evaporator with fan, a thermostatic expansion valve and a coaxial coil heat exchanger as condenser.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003