

Product Name :
Bench Top Mechanical Heat Pump Trainer

Product Code :
NLAB-MECHANICALAB170003



Description :

Bench Top Mechanical Heat Pump Trainer

Technical Specification :

Bench Top Mechanical Heat Pump Trainer: refrigeration and air conditioning lab equipment: air-to-water heat pump trainer consists of a compressor, an evaporator with fan, a thermostatic expansion valve and a coaxial coil heat exchanger as condenser.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003