

**Product Name :**  
Benchtop Fumehood

**Product Code :**  
NLAB-SCIENTIFICEQ12006



**Description :**

Benchtop Fumehood

**Technical Specification :**

Style: Benchtop Ductless

Internal Dimensions (mm): 1600x800x800

Internal volume of the cabinet: 1.0 m3

Display: LCD display: that shows the following information:

Air speed continuously monitored.

Type of filter installed.

Working hours.

Expiration date

Next revision date.

Open door warning

Countdown timer.

Clock and calendar

Filter Columns: 2

Filter number: 2

Number of fans: 2

Average volume of treated air/hour: 160 m<sup>3</sup>/h

Average face velocity: between 0.50 m/s

Conformance to: ASTM E84, NFPA45, ASHARE 110, ANS1295

Lighting: LED

LED intensity: 900 Lux / 2x LED 16 W

Alarms and Scheduled warnings:

Open door warning

Open door in off mode warning

Filter use warning

Next validation warning

Few hours of filter life warning

Countdown timer warning

Expired filter alarm (by hours and date)

Temperature alarm

Equipment without filter alarm.

Low barrier alarm

Power Connection: 240 V, 50/60 Hz

Provide the warranty period of 24 months

Working space to support Sample Preparation & Analysis.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments