

Product Name :
Benchtop Horizontal Laminar Flow Hood

Product Code :
NLAB-SCIENTIFICEQ12011



Description :

Benchtop Horizontal Laminar Flow Hood

Technical Specification :

Principle: Horizontal laminar flow.

Internal volume: 0.30 m³

Average face velocity: 0.02 to 0.40 m/s

Structure: Cold rolled steel coated with epoxy

Work Surface: White secured glass.

Product category: Biohazard class II A

Filters and Efficiency: HEPA Filter with a flow equalizer.

HEPA filters with electronic chip. Motor fan: Double aspiration centrifugal fans.

Control Panel: Soft touch keyboard with a switch for the LED lamps, UV ultraviolet lights and maintenance system.

Lighting: LED lamps with a prismatic diffusing screen 700 to 900 Lux 1 to 2 x 15-16 W.

LCD Display: displays type of filter installed, working hours, expiration date and next revision date, countdown timer, clock and calendar.

Noise level: <30dBA-<50dBA

Alarms and scheduled warnings.

Power Connection: 240 V, 50/60 Hz

Provide the warranty period of 24 months

Working space to support Sample Preparation & Analysis.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003