

Product Name :
 Bend Test Equipment

Product Code :
 NLAB-VOCATIONALTVET13026



Description :

Bend Test Equipment

Technical Specification :

Capacity

200kN	
100mm	100mm cylinder travel
100mm	100mm rod specimen diameter
100mm	100mm plate specimen thickness

<p>50mm specimen width</p>		
<p>180° rotation</p>		
<p>5, 7.5, 10, 12.5, 15, 17.5, 20, 22.5, 25, 27.5, 30, 32.5, 35, 37.5, 4.0, 4.25, 4.5, 4.75, 5.0, 5.25,</p>		
	<p>5.5, 5.75, 6, 8, 10, 12, 14, 16, 18, 20, 24, 28, 30, 32, 36,</p>	
	<p>40, 42, 48, 50, 54, 56, 60, 64, 66, 70, 72, 75, 80, 84, 88,</p>	
	<p>90, 98, 100, 110, 112, 120, 126, 128, 140, 154, 160, 168,</p>	
	<p>175, 192(mm) and etc. 55 types in all. Customized.</p>	
<p>2.2kW/AC 220V/50Hz</p>		
<p>500mm x 860mm x 970mm</p>		
<p>Approx. 1000kg</p>		

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003