


**Product Name :**  
Bend Test Machine

**Product Code :**  
NLAB-MACHINELABWK34016



**Description :**

Bend Test Machine

**Technical Specification :**

To do bending test for re-in forced bar, steel and pipes.

Must be able to do bending testing in both directions according to standard requirements

Must conform to ADTM and ISO requirements.

Type: Motorised

Diameter range of bars and pipes: 50mm

Bending angle – Positive direction: 0o to 180o

Bending angle – Reverse direction: 0o to 25o

Rotating Speed: <37 rpm

Diameter of driving tray: 580mm

Power of Motor: 1.5kW

Power Source: 220V 50Hz, single phase

Must come complete with all standard accessories and spares (rollers, mandrels, etc.)

Standard Welding Equipment Laboratory Set Bend Test Machine Manufacturer, Supplier and Bend Test Machine Exporter in India

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003