

**Product Name :**  
Bending Moment Apparatus

**Product Code :**  
NLAB-MECHANICALAB240001



**Description :**

Bending Moment Apparatus

**Technical Specification :**

Bending Moment Apparatus: Applied Mechanics lab equipment: For measurement of bending moment at a section in a beam. The beam has two parts with a ball bearing hinge at the joint. An under hung spring balance with an adjustable length holds the hinge in compression, hence resists bending moment. The beam is simply supported on two end bearings, one with an adjustable height.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003