


Product Name :
Bending Moment Apparatus

Product Code :
NLAB-MECHANICALAB240001



Description :

Bending Moment Apparatus

Technical Specification :

Bending Moment Apparatus: Applied Mechanics lab equipment: For measurement of bending moment at a section in a beam. The beam has two parts with a ball bearing hinge at the joint. An under hung spring balance with an adjustable length holds the hinge in compression, hence resists bending moment. The beam is simply supported on two end bearings, one with an adjustable height.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003