

Product Name :
Bending of Beam Apparatus

Product Code :
NLAB-MECHANICALAB240033



Description :

Bending of Beam Apparatus

Technical Specification :

Consisting of two G.clamp with knife-edge top. Comprising a metal bar of 100 cm length having width 25 mm and thickness 4.5 mm. Sharp knife edge with hanger is also provided which moves on metal bar. A spherometer fitted on a cast iron stand (Specially designed for this apparatus) with electric contacts is also provided, but without weight.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003