

Product Name :
Nauga Lab Benedict Roth Recording Spirometer

Product Code :
NLAB-LABORATORYEM21012



Description :

Nauga Lab Benedict Roth Recording Spirometer

Technical Specification :

Benedict Roth Recording Spirometer (BMR) consists of following features:-

Six
litre capacity, Spirometer has been provided with a unit of 4-speed
electric recording Kymograph with gravity writing ink-pen.

Easily accessible valves, soda-lime container with screw connections in central chamber drain cocks to
all tubes and containers.

Sampling cock provided for connecting the patient to Spirometer or atmosphere

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments