

**Product Name :**  
Benkelman Beam

**Product Code :**  
NLAB-BTMENLAB150014



**Description :**

Benkelman Beam

**Technical Specification :**

Lightweight Aluminium construction

Ease of Transportation

Unique Telescopic Design Simplifying Field set up

Compact, Thereby reducing the amount of storage space needed

Benkelman Beam

Benkelman Beam utilizes the technique of using balanced beam in conjunction with a suitable vehicle to measure road flexure. The improved Benkelman Beam is a convenient, accurate device for measuring the deflection of flexible pavements under moving wheel loads. Operating on a simple lever arm principle, the unit consists. Supplied with carrying case

Note: Benkelman Beam with Digital Dial Gauge also available at an extra cost

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003