

**Product Name :**  
Bennewart Flex Tester Machines

**Product Code :**  
NLAB-TESTINGLABM83001



**Description :**

Bennewart Flex Tester Machines

**Technical Specification :**

Material: Stainless Steel

Electric Source: 220V

Accuracy: 1%

Weight: 22kg

Frequency: 50/60 Hz

Flexing Speed: 15 in/s

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003