

Product Name :
Bernoulli Tubes

Product Code :
NLAB-EDUSCHOOLEQ152002



Description :

Bernoulli Tubes

Technical Specification :

Three side arms extend from the horizontal tube and allow the attachment of glass tubes to form simple manometers.

For demonstrating that when a fluid passes through a tube and reaches a constriction the fluid speed increases with a loss in pressure.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003