

Product Name :
Bicycle Ergograph

Product Code :
NLAB-EDUSCHOOLEQ28001



Description :

Bicycle Ergograph

Technical Specification :

Additionally,
customers are eased with the availability of the Wooden Stand Bicycle
Ergograph loaded with all the advanced features.

A bicycle frame is supported by a wooden stand, from the front of which two uprights ascend and carry a desk and cross piece.

Which
provide for the attachment, directly in front of the subject, of the
tension balances and other pieces in front of the subject, of the
tension balances and other pieces of apparatus.

A cast iron wheel of 5 1/2 feet in circumference, weighing 22 Kilos and mounted on

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003