


Product Name :
Bifilar And Trifilar Apparatus

Product Code :
NLAB-MECHANICALAB240003



Description :

Bifilar And Trifilar Apparatus

Technical Specification :

Bifilar And Trifilar Apparatus: In a bifilar suspension, the pendulum body is suspended over two threads. The pendulum body oscillates in a plane purely transnationally without rotation. This kind of pendulum can be considered as a mathematical pendulum. In a trifiler suspension with three threads, the pendulum body is set in a torsional vibration. The torsional vibration can be used

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003