

Product Name : Bile Broth Base	Product Code : NLAB-CHEMICALSIND21027
Description :	
Bile Broth Base	
Technical Specification :	
Bile Broth Base	
Bile Broth Base is used for the cultivation of members of Enterobacteriaceae and in culture of blood clots from patients with suspected enteric fever.	
Directions	
Suspend 30.0 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 40°C and add 1 ml of Streptokinase solution. Mix well and dispense as desired.	
Quality Control	
Reaction: Reaction of 3% w/v aqueous solution at 25°C. pH : 7.6±0.2	
pH: 7.40-7.80	
Appearance: Cream to yellow homogeneous free flowing powder	
Colour and Clarity of prepared medium: Yellow coloured, clear solution without any haziness.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003