

Product Name :
Binocular Microscope

Product Code :
NLAB-ELECTRONICSAB150002



Description :

Binocular Microscope

Technical Specification :

Binocular Microscope Science Lab
Head rotated 360 degrees.
Transmitted illumination.
Adjustable interpupillary distance between 55mm- 75mm.
Includes Spare Bulb, Fuse, Power Cord and Plastic Dust Cover.
Ideal for general medical applications, industrial use, hobbyists/ collectors and educational purposes.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003