

Product Name :
Binocular Microscope

Product Code :
NLAB-EDUSCHOOLEQ173076

Description :

Binocular Microscope

Technical Specification :

The presence of co-axial coarse as well as fine focusing adjustment knobs also provides for fine motion that can be graduated to .002 mm.

Binocular Microscope offered come in 45° inclined finish with 360 degree rotatable binocular tube.

Further, the steeple grey as well as chemical resistant back-on finish and heavy rectangular sturdy base that comes with built-in illumination through 6 V - 20 W halogen lamp also provides for easy usage.

For smooth movement, it also features presence of dust proof quadruple/quintuple ball bearing for meeting the revolving needs of nose piece with positive click stops.

Features:-

Exclusive preset focus lock prevents damage to valuable slides and objectives.

Heavy rectangular sturdy base, with built-in illumination 6 V - 20 W halogen lamp.

45° inclined with 360 degree rotatable binocular tube.

Mechanical Tube Length 160 mm.

Dust proof quadruple/quintuple ball bearing revolving nose piece with positive click stops.

Steeple Grey, chemical resistant, back-on finish.

Co-axial built in Mechanical stage with fine vernier graduation, designed with convenient coaxial adjustment for Jerk free smooth slide manipulation through 50 X 75 mm.

Co-axial coarse & fine focusing adjustment knobs. Fine motion, graduated to .002 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003