

**Product Name :**  
Binocular Research Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ1730143

**Description :**

Binocular Research Microscope

**Technical Specification :**

Nosepiece: Quadruple nosepiece with accurate centering and positive click stops.

Eyepieces: Wide field WF 10X and H5x or H6X or H15X (any two pairs).

Observation Head: 45° inclined Binocular head rotatable through 360°.

Stage: Built in graduated mechanical stage of 135x120mm. Low coaxial drive controls for easy and smooth scan of specimen slide over range of 50x75mm.

Illumination: 6V-20W halogen lamp, 220/110Volts and mirror attachment for working in day light.

Objectives: Achromatic 4X, 10X, 40X SL, 100X SL Oil immersion.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003