

**Product Name :**  
Binocular Stereo Zoom Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ1730138

**Description :**

Binocular Stereo Zoom Microscope

**Technical Specification :**

Pillar style table stand has reversible black and white stage plate to provide contrast with light- and dark-colored specimens, and large stage clips to secure the specimen during viewing.

Professional stereo microscope with pillar stand has long working distance to enable users to perform work or manipulate large items, including circuit boards and dental appliances.

Binocular viewing head with pair of 10x super-widefield eyepieces, adjustable interpupillary distance, fixed 45-degree vertical inclination to reduce eye and neck strain, and 360-degree rotation capability to enable sharing.

Removable rheostat-controlled 64-bulb LED ring light provides intensity control and bright, cool light for working with temperature-sensitive or live specimens.

0.7x-4.5x zoom objective provides continuous zoom magnification and longer focal length for inspecting large-scale specimens, and a 0.5x Barlow lens extends the working distance.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003