

Product Name : Binocular Upright Metallurgical Microscope	Product Code : NLAB-EDUSCHOOLEQ1730142
Description :	
Binocular Upright Metallurgical Microscope	
Technical Specification :	
Nose piece: Quadruple revolving nose piece with accurate centering & positive click stops.	
Eyepieces: H6X-WF10X-H16X (any two).	
Body: Binocular Head inclined at 45-degrees.	
Objectives: Achromatic DIN M4X-M10X & M40X.	
Focusing: Both side co-axial focusing knobs.	
Illumination: Brightfield vertical illuminator 6V-20W Halogen with field diaphragm, filter slots and bulb centring mechanism.	
Stage: Large Co-axial Mechanical Stage for X-Y movement of specimens.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003