


**Product Name :**  
Biology Lab Clinostat Electric

**Product Code :**  
NLAB-EDUSCHOOLEQ13005



**Description :**

Biology Lab Clinostat Electric

**Technical Specification :**

The Clinostat is supplied with a 70mm diameter cork lined/faced disc with a fitting to take a specimen tube.

Essential for studies on geotropism and phototropism on seedlings, cut shoots, and small potted plants.

The angle is adjustable by a thumb screw with a clear plastic cover.

The motor rotates 5 revolutions per hour.

220V, 50Hz. 300mm (H).

The seedlings are pinned to the cork table which is rotated by an electric motor.

Disc diameter: 100mm.

Complete with 3 core cable.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003