

Product Name :
Biology Lab Concentrated Enzyme Buffer Ph 2.8

Product Code :
NLAB-EDUSCHOOLEQ22006



Description :

Biology Lab Concentrated Enzyme Buffer Ph 2.8

Technical Specification :

Enzyme buffer solutions are presented as 50ml concentrates in handy storage bottles.

These buffers are meticulously prepared by experienced technicians working in clean laboratory production facilities.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003