

Product Name :
Biology Lab Ganong Potometer Apparatus

Product Code :
NLAB-EDUSCHOOLEQ13032



Description :

Biology Lab Ganong Potometer Apparatus

Technical Specification :

A plant shoot is sealed into a stopper which fits the end tube and the whole apparatus is filled with water.

For the demonstration of transpiration rates and rates of water absorption by a cut shoot under various external conditions.

The position of the air bubble may easily be adjusted by opening the glass stopcock to admit water from the reservoir.

The open end of the capillary tube is dipped into a beaker of water.

The whole apparatus is mounted on a varnished wooden base by spring clips.

Capillary Bore: 1mm.

Base: 20 x 80mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003