

Product Name :
Biology Lab Poppit Beads Red

Product Code :
NLAB-EDUSCHOOLEQ23008



Description :

Biology Lab Poppit Beads Red

Technical Specification :

They can represent disaccharide and polysaccharide molecules in demonstrating carbohydrate metabolism.

Poppit beads are a useful laboratory modelling unit.

Poppit Beads can also be used to represent genetic code as well as whole chromosomes, or as dipeptides, amino acids, monosaccharide, polypeptides and proteins in protein metabolism.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003