

Product Name :
Biology Lab Poppit Beads Yellow

Product Code :
NLAB-EDUSCHOOLEQ23009



Description :

Biology Lab Poppit Beads Yellow

Technical Specification :

Poppit beads are a useful laboratory modelling unit.

They can represent disaccharide and polysaccharide molecules in demonstrating carbohydrate metabolism.

Poppit Beads can also be used to represent genetic code as well as whole chromosomes, or as dipeptides, amino acids, monosaccharide, polypeptides and proteins in protein metabolism.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003