

Product Name :
Biology Lab Potometer Apparatus H Type

Product Code :
NLAB-EDUSCHOOLEQ13031



Description :

Biology Lab Potometer Apparatus H Type

Technical Specification :

The cut shoot is attached to the side arm with tubing and the wide glass tube acts as a reservoir.

Suitable for individual student use.

Capillary tube length 140mm.

Glass tube length 160 x 20mm Dia.

Bore: 1mm.

Side arm 8mm dia.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

