


Product Name :
Biology Lab Potometer Apparatus White Stand

Product Code :
NLAB-EDUSCHOOLEQ13030



Description :

Biology Lab Potometer Apparatus White Stand

Technical Specification :

A syringe is used to zero the scale.

When transpiration occurs, the bubble within the capillary moves.

The apparatus is submerged and the specimen placed in the tube.

Used to show transpiration which occurs in small plant shoots.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003