


**Product Name :**  
Biology Lab Tulgren Funnel

**Product Code :**  
NLAB-EDUSCHOOLEQ16012



**Description :**

Biology Lab Tulgren Funnel

**Technical Specification :**

A soil sample is placed on the perforated plate.

A 220V heating lamp is held above the soil, driving the organisms down through the funnel into a collecting jar.

Supplied with coarse and fine mesh and 3 core cable.

Requires a low wattage bulb and a beaker (not provided).

Retort Rod and Base not included.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003