

Product Name :
Bitumen and Centrifuge Extractor - Motorized

Product Code :
NLAB-BTMENLAB150002



Description :

Bitumen and Centrifuge Extractor - Motorized

Technical Specification :

Same as above but shaft is rotated by an electric motor and gear. Care is taken to prevent solvent entering into the rotor of electrical motor. Operates on 230 Volts A.C. single phase.

Accessories : 1. Filter paper disc (packet of 100nos.) 2. Benzene

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003