

**Product Name :**  
BJT Biasing Trainer Digital Meters

**Product Code :**  
NLAB-ELECTRICALAB210032



**Description :**

BJT Biasing Trainer Digital Meters

**Technical Specification :**

**SCOPE OF LEARNING:**

Fixed Biasing Without Emitter Resistor

Collector To Base Bias

Fixed Biasing With Emitter Resistor

Potential Divider Biasing/ Self Biasing

**TECHNICAL SPECIFICATIONS:**

Digital Meters:

Voltmeter 20V/2V DC.

Ammeter 200mA/20mA DC.

Ammeter 200uA/2mA DC.

Power Supplies:

DC Supply IC Regulated +12V DC, 250mA.

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

Transistor SL100

Required Resistors

SALIENT FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Connections are brought out through 2mm Colored Sockets.

Patch Cords 2mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x10"

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003