

Product Name :
Blaine Permeability Meter for Fineness Determination

Product Code :
NLAB-EDUCATIONTVET15042



Description :

Blaine Permeability Meter for Fineness Determination

Technical Specification :

According to EN 196/6, ASTM C204 - AASHTOT153, BS4359 / 2, UNI 7374, NF P15: 442, UNE80106

Reference cement

Manometric liquid 200 ml

Filter paper (1000 pcs)

Reference sand I (large)

Reference sand II (fine)

0.01 mm digital caliper

Chattaway spatula

Precision balance, type PR 423 / E

Basic weighing (10 units), piece counting, weighing in%

Backlit liquid crystal display (LCD)

Communication interface: RS-232

Spirit level, leveling feet

Weighing chamber with 3 transparent sliding doors

Weighing hook under scale

Stainless steel weighing pan: Ø 120 mm

External calibration

Weight: 420 g

Display resolution: 1 mg

Repeatability: 1 mg

Power supply: 230V - 50Hz

Calibration mass, class F1, 200 g, delivered with certificate

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003