

Product Name :
Blast Freezer / Chest Freezer

Product Code :
NLAB-SCIENTIFICEQ35004



Description :

Blast Freezer / Chest Freezer

Technical Specification :

Blast Freezer are used for instant freezing of blood plasma, food products etc. it freezes the blood plasma and other samples to -20°C in not less than 60 min. Chest Freezer offer a breakthrough in sample preparation for use prior to cold storage including flexible design, accelerated air flow and a self-contained refrigeration system. Slow freezing forms big ice crystals in the product. The crystals are forming cells and destroy the structure of the product but in blast freezer Fast freezing forms many small ice crystals which keeps the original product safe.

- 93/42/EEC directive and amendments
- BS 4376-2:1996

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003