

**Product Name :**  
Blast Freezer / Chest Freezer

**Product Code :**  
NLAB-SCIENTIFICEQ35004



**Description :**

Blast Freezer / Chest Freezer

**Technical Specification :**

Blast Freezer are used for instant freezing of blood plasma, food products etc. it freezes the blood plasma and other samples to -20oC in not less than 60 min. Chest Freezer offer a breakthrough insample preparation for use prior to cold storage including flexible design, accelerated air flow and a self-contained refrigeration system. Slow freezing forms big ice crystals in the product. The crystals are forming cells and destroy the structure of the product but in blast freezer Fast freezing forms many small ice crystals which keeps the original product safe.

- 93/42/EEC directive and amendments
- BS 4376-2:1996

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003