

**Product Name :**  
Bleeding Of Freshly Mixed Concrete

**Product Code :**  
NLAB -CONCRLAB120035



**Description :**

Bleeding Of Freshly Mixed Concrete

**Technical Specification :**

Stainless steel cylindrical container for bleeding test to EN 480-4, approx. 14 liters capacity.

Internal dia 255 +/- 5mm x 280 +/- 5mm inside height. Complete with stainless steel cover. Flat rounded steel float with a diameter of  $(100 \pm 10)$  mm.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003