

Product Name :
Blood Bag Double CPDA-1

Product Code :
NLABMEDICAL-M2290001



Description :

Blood Bag Double CPDA-1

Technical Specification :

Blood Bag Double CPDA-1

Basic Configuration:

- * Blow Extruded
- * Rolled
- * Donor Needle (16G/17G/18G/20G)
- * Platelet storage bag (for platelet storage 5 or 7 days)
- * Sampling Device:- Needle Guard
- * Sampling Device:- Vacuum Tube Holder
- * Sampling Device:- Vacuum tube holder and sampling bag

We are Exporters and Suppliers of Blood Bag Double CPDA-1 at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003