

|  |   |
|--|---|
| <b>Product Name :</b><br>Blood Bag Sealer  | <b>Product Code :</b><br>NLABMEDICAL-M3200008 |
|  |   |
| <b>Description :</b>   |   |
| Blood Bag Sealer India, Blood Bag Sealer manufacturer India, suppliers, Medical Blood Bag Sealer India, Hospital Blood Bag Sealer India    |   |
| <b>Technical Specification :</b>   |   |
| Sealing operation is activated automatically when the tube is based in between the electrodes.   |   |
| Electrodes are made from high grade metal for perfect sealing. Caps provided in front of electrodes are removable, so it is easy to clean. |   |
| LED's are provided for indications of Power, Read, Seal, Ground.   |   |

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003